



530 Main Street, Acton, MA 01720

Phone: (978)263-3584, **Fax:** (978)263-5086

Web Site: www.acton-research.com

**Instructions
for
MonoControl SDK**

SpectraPro OLE COM Interface

The SpectraPro COM interface allows users to create software applications that interface with and control the software and hardware of the Acton Research Corporation SpectraPro monochromators through Microsoft's COM/DCOM protocol. This interface is declared with GUID {999DF143-A7CF-4815-9BBF-D2E71E0929E5} and recognized by the system as *ARCSpectraPro.SpectraPro_Interface*. The interface is self registering on the invocation of SpectraPro.

The Interface is divided into two parts:

1. Program Control : Functions that control/modify how the software appears on the screen.
2. Hardware Control : Functions that control/modify the Acton hardware in the system.

A complete set of example programs written in Visual Basic, and Delphi are included on the SpectraPro installation disk.

Program Control Functions

1. SetPageNum Brings the specified SpectraPro page to the front
2. GetPageNum Returns the current SpectraPro screen
3. IsSpectraProUp Returns SpectraPro loading status
4. SpectraProToFront Brings SpectraPro to top of desktop
5. SpectraProToBack Sends SpectraPro to bottom position on desktop
6. SpectraProShow Restores SpectraPro to the Screen
7. SpectraProHide Removes SpectraPro from the Screen
8. getSpectraProMono When used with the NCL Data Acquisition Module, returns the number of the mono with which the COM interface is communicating (1 or 2).
9. setSpectraProMono When used with the NCL Data Acquisition Module, sets the mono (1 or 2) with which the SpectraPro COM interface is communicating.
10. UpdateSpectraPro Refreshes the SpectraPro screen with the latest information

SetPageNum

Description: This function will select which of the SpectraPro pages is to be visible.

Parameters

PageNum:	0	About Screen
	1	Main Page
	2	Operation Page
	3	Diverter Page
	4	Motorized Slits Page
	5	Install Gratings Page
	6	Wave Offset Page
	7	G Adjust Page
	8	Monochromator Defaults Page
	9	Calibration Info Page
	10	Terminal Page

Result:	0	An invalid request
	!0	A valid request

Definitions

Delphi

Function SetPageNum(PageNum: Integer): Integer; safecall;

Microsoft IDL

```
[id(0x00000007)]  
HRESULT SetPageNum(  
    [in] long Page,  
    [out, retval] long* ValidRequest);
```

Visual Basic Sample Code

```
Dim PageNum As Long  
  
If SpectraPro_COM.SelectPage(PageNum) = 0 then  
Msgbox ("Failed to select page")  
End If
```

Delphi Sample Code

```
Var  
PageNum : Long;  
  
Begin  
SpectraProServer.SelectPage(PageNum) = 0;  
Then show message ('Failed to select page');  
End;
```

getPageNum

Description: This function will return which of the SpectraPro pages is currently up.

Parameters

Result:

- 0 About Screen
- 1 Main Page
- 2 Operation Page
- 3 Diverter Page
- 4 Motorized Slits Page
- 5 Install Gratings Page
- 6 Wave Offset Page
- 7 G Adjust Page
- 8 Monochromator Defaults Page
- 9 Calibration Info Page
- 10 Terminal Page

Definitions

Delphi

Function GetPageNum: Integer; safecall;

Microsoft IDL

[id(0x00000008)]
HRESULT getPageNum([out, retval] long* Page);

Visual Basic Sample Code

```
Dim PageNum As Long  
PageNum = SpectraPro_COM.getPageNum
```

Delphi Sample Code

```
Var  
PageNum : Long;  
  
Begin  
PageNum := SpectraProServer.getPageNum;  
End;
```

IsSpectraProUp

Definition: Tells if SpectraPro has finished loading

Parameters

Result: 0 has not finished
!0 has finished

Definitions

Delphi

Function IsSpectraProUp: Integer; safecall;

Microsoft IDL

[id(0x00000001)]
HRESULT IsSpectraProUp([out, retval] long* SpectraProState);

Visual Basic Sample Code

```
If SpectraPro_COM.IsSpectraProUp = 0 then  
    Labell.caption="SpectraPro is loading"  
End if
```

Delphi Example Code

```
if SpectraProServer.IsSpectraProUp = 0  
    then label.caption = 'SpectraPro is loading';
```

SpectraProToFront

Description: Brings SpectraPro to the top of the desktop

Definitions

Delphi

Procedure SpectraProToFront; safecall;

Microsoft IDL

[id(0x00000006)]
HRESULT SpectraProToFront();

Visual Basic Sample Code

```
SpectraPro_COM.SpectraProToFront
```

Delphi Sample Code

```
SpectraProServer.SpectraProToFront;
```

SpectraProToBack

Description: Sends SpectraPro to the bottom of the desktop

Definitions

Delphi

Procedure SpectraProToBack; safecall;

Microsoft IDL

[id(0x00000005)]
HRESULT SpectraProToBack();

Visual Basic Sample Code

```
SpectraPro_COM.SpectraProToBack
```

Delphi Sample Code

```
SpectraProServer.SpectraProToBack;
```


SpectraProShow

Description: Display's SpectraPro on the Screen

Definitions

Delphi

Procedure SpectraProShow; safecall;

Microsoft IDL

[id(0x00000002)]
HRESULT SpectraProShow();

Visual Basic Sample Code

```
SpectraPro_COM.SpectraProShow
```

Delphi Sample Code

```
SpectraProServer.SpectraProShow;
```

SpectraProHide

Description: Remove SpectraPro from the screen

Definitions

Delphi

Procedure SpectraProShow; safecall;

Microsoft IDL

[id(0x00000003)]
HRESULT SpectraProHide();

Visual Basic Sample Code

```
SpectraPro_COM.SpectraProHide
```

Delphi Sample Code

```
SpectraProServer.SpectraProHide;
```

getSpectraProMono

Description: retrieves the number of the monochromator (1 or 2) with which the SpectraPro COM interface software is communicating.

Parameters

Result: The number of the monochromator with which the SpectraPro COM interface is communicating can be 1 or 2 if an NCL is present. Otherwise a 1 is returned.

Definitions

Delphi

```
function getSpectraProMono: Integer; safecall;
```

Microsoft IDL

```
[id(0x0000001b)]  
HRESULT getSpectraProMono([out, retval] long* curMono);
```

Visual Basic Sample Code

```
Dim MonoNum as Long  
MonoNum = SpectraPro_COM.getSpectraProMono
```

Delphi Sample Code

```
Var MonoNum : integer;  
MonoNum := SpectraProServer.getSpectraProMono;
```

setSpectraProMono

Description: Sets which Monochromator with which the SpectraPro COM interface software is communicating.

Parameters

NewMono: The number of the monochromator to be selected, 1 valid for singles mono's, 1 and 2 valid for mono's with NCL's.

Result: 0, Operation failed.
!0, Operation succeeded.

Definitions

Delphi

```
function setSpectraProMono(newMono: Integer): Integer; safecall;
```

Microsoft IDL

```
[id(0x0000001c)]  
HRESULT setSpectraProMono(  
    [in] long newMono,  
    [out, retval] long* ValidResult);
```

Visual Basic Sample Code

```
Dim MonoNum as Long  
  
If SpectraPro_COM.setSpectraProMono(MonoNum) = 0 then  
    MsgBox ("Failed")  
End if
```

Delphi Sample Code

```
Var MonoNum : integer;  
  
If SpectraProServer.setSpectraProMono(MonoNum) <> 0  
    Then showmessage('Failed');
```

UpdateSpectraPro

This function is used to update the values on the SpectraPro software screens after changes have been initiated in your application. This function is especially useful in debugging your code as it will verify if you have correctly programmed the instrument state, or input values.

Description: Low level changes to SpectraPro are not always immediately translated to the SpectraPro screens. Run this function to resynchronize the screens to the current conditions.

Definitions

Delphi

```
Procedure UpdateSpectraPro; safecall;
```

Microsoft IDL

```
[id(0x00000004)]  
HRESULT UpdateSpectraPro();
```

Visual Basic Sample Code

```
SpectraPro_COM.UpdateSpectraPro
```

Delphi Sample Code

```
SpectraProServer.UpdateSpectraPro;
```

SpectraPro OLE COM interface Hardware Control Functions

Monochromators

1. Mono_Present	Tests to see if a monochromator is in the system
2. MonoDouble	Defines a monochromator as a double monochromator
3. getMonoWavelength	Returns a monochromator position in nm
4. setMonoWavelength	Sends the monochromator to a wavelength in nm
5. getMonoGrating	Returns the grating number in use
6. setMonoGrating	Change the grating to the number specified
7. getMonoTurret	Returns the turret number in use
8. setMonoTurret	Change the turret to the number specified
9. getMonoDiverterPosition	Returns the position of the active port
10. setMonoDiverterPosition	Change the diverter mirror position
11. MonoSlitMotorized	Tests to see if a slit is motorized or manual
12. getMonoSlitWidth	Returns the slit width of a motorized slit
13. setMonoSlitWidth	set a motorized slit to a specific width
14. MonoSetupScan	Setup the parameters for scanning a Monochromator
15. getMonoScanParams	Find out the scan parameters for a Monochromator
16. MonoScanWavelength	Scan a Monochromator

NCL

- 1. NCL_Present**

Filters

1. Filter_present	Test to see if a filter wheel is present
2. getFilterPosition	Returns the current filter position (1 – 6)
3. setFilterPosition	Moves the specified filter into position
4. FilterHome	Homes the filter wheel and sets it to position 1

Mono_Present

Description: Used to determine if a specific monochromator is in the system. When used with an NCL, use the parameters 1 or 2 to determine if one or both monochromators are present in the system. When used without an NCL, use input parameter 1 to determine if the monochromator is actively attached to the system.

Parameters

Mono_Num The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

Result:

0	Not Present
!0	Present

Definitions

Delphi

```
Function Mono_Present( Mono_Num:Integer): Integer; safecall;
```

Microsoft IDL

```
[id(0x0000000a)]  
HRESULT Mono_Present(  
    [in] long MonoNum,  
    [out, retval] long* MonoState);
```

Visual Basic Example Code

```
If SpectraPro_COM.Mono_Present(1) = 0 then  
    Labell.caption = "Monochromator 1 not present"  
End if
```

Delphi Example Code

```
if SpectraProServer.Mon_Present(2) = 0  
    then labell.caption = 'Mono 2 not present';
```

MonoDouble

Description: Determines if a monochromator is a double or single

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

Result: 0 if it is not a double monochromator
!0 if it is a double monochromator

Definitions

Delphi

```
Function MonoDouble( Mono_Num: Integer  
): Integer; safecall;
```

Microsoft IDL

```
[id(0x0000000c)]  
HRESULT MonoDouble(  
    [in] long MonoNum,  
    [out, retval] long* DoubleState);
```

Visual Basic Sample Code

```
If SpectraPro_COM.MonoDouble(1) = 0 then  
    Labell.caption = "Monochromator 1 is not a double"  
End if
```

Delphi Sample Code

```
if SpectraProServer.MonDouble(2) <> 0  
    then Showmessage ('Monochromator 2 is a double');
```


getMonoWavelength

Description: This function will return the position of the specified monochromator in nm

Parameters

Mono_Num: The number of the monochromator being interrogated (see note below).

Result: The monochromator position in nm

Definitions

Delphi

```
Function getMonoWavelength( Mono_Num:Integer  
                           ): Double; safecall;
```

Microsoft IDL

```
[id(0x0000000d)]  
HRESULT getMonoWaveLength(  
    [in] long MonoNum,  
    [out, retval] double* CurWave);
```

Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

Visual Basic Sample Code

```
Dim Monowavelength As Double  
  
Monowavelength = SpectraPro_COM.getMonoWavelength(1)  
Label1.Caption = "Mono 1 is at" & Monowavelength & " nm"
```

Delphi Sample Code

```
Var MonoWavelength : double  
  
MonoWavelength := SpectraProServer.getMonoWavelength(1);  
Label1.caption := 'Wavelength : ' +format('% .3f',[TempDbl]) + ' nm';
```

setMonoWaveLength

Description: This function will change the position of the specified monochromator.

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the connection between the monochromator and NCL (1 or 2).

NewWavelength: Enter the new position for the specified monochromator in NM.
Note: NewWavelength may be integer or decimal values.

Result:
0 An invalid request
!0 A valid request

Definitions

Delphi

```
function setMonoWaveLength(      Mono_Num:Integer;
                               NewWaveLength: Double
                               ): Integer; safecall
```

Microsoft IDL

```
[id(0x0000000e)]
HRESULT setMonoWaveLength(
    [in] long MonoNum,
    [in] double newWave,
    [out, retval] long* ValidRequest);
```

Visual Basic Sample Code

```
Dim newWave As Double
Dim curmono As Double
' set the wavelength
Curmono = 1
newWave = 500.0
If SpectraPro_COM.setMonoWavelength(curmono, newWave) = 0 Then
    MsgBox ("Error : Failed to set wavelength")
End If
```

Delphi Example Code

```
var newWave : double;
var curmon : double;
// set the wavelength
newWave := 500.0;
curmono := 1;

if SpectraProServer.setMonoWavelength(curmono,newWave) = 0;
  then ShowMessage('Error : Failed to set wavelength');
```

getMonoGrating

Description: Returns the currently selected grating number

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

Result: The current grating number in the specified monochromator

Definitions

Delphi

```
Function getMonoGrating( Mono_Num:integer;  
                        ): Integer; safecall;
```

Microsoft IDL

```
[id(0x0000000f)]  
HRESULT getMonoGrating(  
    [in] long MonoNum,  
    [out, retval] long* curGrating);
```

.

Visual Basic Sample Code

```
Dim newGrat As Integer  
  
newGrat = SpectraPro_COM.getMonoGrating(1)  
Labell.caption = newGrat
```

Delphi Sample Code

```
var GratingNumber : integer;  
  
GratingNumber := SpectraProServer.getMonoGrating(1)  
ShowMessage('Mono 1 grating is : ' +inttoStr(grating number));
```

setMonoGrating

Description: Puts the specified grating number in position for the specified monochromator

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

New Grating: The grating number to be put in place

Result: 0 An invalid request
!0 A valid request

Definitions

Delphi

```
Function setMonoGrating( Mono_Num: Integer;  
                        NewGrating: Integer  
                        ): Integer; safecall;
```

Microsoft IDL

```
[id(0x00000010)]  
HRESULT setMonoGrating(  
    [in] long MonoNum,  
    [in] long newGrating,  
    [out, retval] long* ValidRequest);
```

Visual Basic Sample Code

```
Dim newGrat As Long  
  
newGrat = 3  
' set the grating  
If SpectraPro_COM.setMonoGrating(1, newGrat) = 0 Then  
    MsgBox ("Error : Failed to move grating")  
End If
```

Delphi Sample Code

```
var newGrat : integer;  
  
newGrat := 2;  
// set the grating  
if SpectraProServer.setMonoGrating(1,newGrat) = 0  
    then ShowMessage('Error : Failed to move grating');
```

getMonoTurret

Description: Returns the currently selected grating turret number

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

Result: The current turret number in the specified monochromator

Definitions

Delphi

```
Function getMonoTurret( Mono_Num:integer;  
                        ): Integer; safecall;
```

Microsoft IDL

```
[id(0x00000011)]  
HRESULT getMonoTurret(  
    [in] long MonoNum,  
    [out, retval] long* curTur);
```

Visual Basic Sample Code

```
Dim newTur As Integer  
  
newTur = SpectraPro_COM.getMonoTurret(1)  
Label1.caption = newTur
```

Delphi Sample Code

```
var TurretNumber : integer;  
  
TurretNumber := SpectraProServer.getMonoTurret(1)  
ShowMessage('Mono 1 turret is : ' +inttoStr(TurretNumber));
```

setMonoTurret

Description: selects the specified turret for the monochromator

Parameters

Mono_Num:
New Grating: The turret number to be selected

Result:
0 An invalid request
!0 A valid request

Definitions

Delphi

```
Function setMonoTurret( Mono_Num: Integer;  
NewTurret: Integer  
): Integer; safecall;
```

Microsoft IDL

```
[id(0x00000012)]  
HRESULT setMonoTurret(  
[in] long MonoNum,  
[in] long newTur,  
[out, retval] long* ValidRequest);
```

Visual Basic Sample Code

```
Dim newTur As Long  
  
newTur = 3  
' set the turret  
If SpectraPro_COM.setMonoGrating(1, newTur) = 0 Then  
MsgBox ("Error : Failed to change turret")  
End If
```

Delphi Sample Code

```
var newTurret : integer;  
  
newTurret := 2;  
// set the grating  
if SpectraProServer.setMonoGrating(1,newTurret) = 0  
then ShowMessage('Error : Failed to change turret');
```

getMonoDiverterPosition

Description: This function indicates which slit is in the light path. Note that the same function is used with different variables for determining the active entrance and exit slit.

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

Divert_Num:

1	Entrance mirror
2	Exit mirror
3	Slave Entry Mirror (double monochromator)
4	Slave Exit mirror (double monochromator)

Result:

0	Failed
1	side entrance slit
2	front entrance slit
3	front exit slit
4	side exit slit
5	slave side entrance slit
6	slave front entrance slit
7	slave front exit slit
8	slave side exit slit

Definitions

Delphi

```
Function getMonoDiverterPosition( Mono_Num: Integer;  
                                Divert_Num: Integer;  
                                ): Integer; safecall;
```

Microsoft IDL

```
[id(0x00000013)]  
HRESULT getMonoDiverterPosition(  
    [in] long MonoNum,  
    [in] long Divert_Num,  
    [out, retval] long* curSlit);
```

Visual Basic Sample Code

```
Dim curmono As Long  
If SpectraPro_COM.setMonoDiverterPosition(curmono, 1) = 1 Then  
    Labell.caption "The side entrance slit is active"  
End If
```


Delphi Sample Code

```
// entrance position, 1 = side, 2 = front  
If SpectraProServer.getMonoDiverterPosition(curmono,1)  
Then showmessage ('Side Entrance')  
else showmessage ('Front Entrance');
```

setMonoDiverterPosition

Description: This function is used to put a specific slit in the light path. This is done by moving the appropriate diverter mirror.

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

Divert_Num:

1	Entrance mirror
2	Exit mirror
1	Slave Entry Mirror (double monochromator)
2	Slave Exit mirror (double monochromator)

NewSlit_Num:

0	Failed
1	side entrance slit
2	front entrance slit
3	front exit slit
4	side exit slit
5	slave side entrance (doubles only usually the intermediate slit)
6	slave front entrance (doubles only)
7	slave front exit slit (doubles only)
8	slave side exit slit (doubles only)

Result:

0	An invalid request
!0	A valid request

Definitions

Delphi

```
Function setMonoDiverterPosition(           Mono_Num: Integer;  
                                       Divert_Num: Integer;  
                                       NewSlit_Num: Integer  
                                       ): Integer; safecall;
```

Microsoft IDL

```
[id(0x00000014)]  
HRESULT setMonoDiverterPosition(  
    [in] long MonoNum,  
    [in] long Divert_Num,  
    [in] long NewSlit,  
    [out, retval] long* ValidRequest);
```

Visual Basic Sample Code

```
If SpectraPro_COM.setMonoDiverterPosition(1, 1, 1) = 0 Then  
    MsgBox ("Error : Failed to change divertor")  
End If
```

Delphi Sample Code

```
var curmono : Integer  
  
// flip entrance diverter  
if SpectraProServer.setMonoDiverterPosition(curmono,1,1) = 0  
    then ShowMessage('Error : Failed to change diverter');
```

MonoSlitMotorized

Description: Determines if a particular slit is manual or motorized

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

Slit_Num:

0	Failed
1	side entrance slit
2	front entrance slit
3	front exit slit
4	side exit slit
5	slave side entrance (doubles only usually the intermediate slit)
6	slave front entrance (doubles only)
7	slave front exit slit (doubles only)
8	slave side exit slit (doubles only)

Result:

0	An invalid request
!0	A valid request

Definitions

Delphi

```
Function MonoSlitMotorized(      Mono_Num: Integer;  
                             Slit_Num: Integer  
                             ): Integer; safecall;
```

Microsoft IDL

```
[id(0x00000015)]  
HRESULT MonoSlitMotorized(  
    [in] long MonoNum,  
    [in] long SlitNum,  
    [out, retval] long* SlitState);
```

Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL.

Visual Basic Sample Code

```
Dim curmono As Long  
If SpectraPro_COM.MonoSlitMotorized(curmono, 1) <> 0 Then  
    MsgBox ("Side Entrance Slit is Motorized")  
End if
```

Delphi Sample Code

```
// side entrance slit , slit # 1
```

```
if SpectraProServer.MonoSlitMotorized(curmono,1) <> 0
  then showmessage ('Side Entrance Slit is Motorized')
  else showmessage ('Side Entrance Slit is not Motorized');
```

getMonoSlitWidth

Description: Returns the slit width in microns only.

Parameter

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

Slit_Num:

0	Failed
1	side entrance slit
2	front entrance slit
3	front exit slit
4	side exit slit
5	slave side entrance (doubles only usually the intermediate slit)
6	slave front entrance (doubles only)
7	slave front exit slit (doubles only)
8	slave side exit slit (doubles only)

Result: The width of the slit in microns

Definitions

Delphi

```
Function getMonoSlitWidth( Mono_Num: Integer;  
                           Slit_Num: Integer  
                           ): Integer; safecall;
```

Microsoft IDL

```
[id(0x00000016)]  
HRESULT getMonoSlitWidth(  
    [in] long MonoNum,  
    [in] long SlitNum,  
    [out, retval] long* SlitWidth);
```

Visual Basic Sample Code

```
Label1.Caption = SpectraPro_COM.getMonoSlitWidth(curmono, 1)
```

Delphi Sample Code

```
Label1.caption := IntToStr(SpectraProServer.getMonoSlitWidth(curmono,1));
```

setMonoSlitWidth

Description: Sets the width of the selected slit in microns. Valid only for motorized slits

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the of the connection between the monochromator and NCL (1 or 2).

Slit_Num:

0	Failed
1	side entrance slit
2	front entrance slit
3	front exit slit
4	side exit slit
5	slave side entrance (doubles only usually the intermediate slit)
6	slave front entrance (doubles only)
7	slave front exit slit (doubles only)
8	slave side exit slit (doubles only)

newWidth: The new slit width in microns

Result:

0	An invalid request
!0	A valid request

Definitions

Delphi

```
Function setMonoSlitWidth(           Mono_Num: Integer;  
                               Slit_Num: Integer;  
                               NewWidth: Integer  
                               ): Integer; safecall;
```

Microsoft IDL

```
[id(0x00000017)]  
HRESULT setMonoSlitWidth(  
    [in] long MonoNum,  
    [in] long SlitNum,  
    [in] long newWidth,  
    [out, retval] long* ValidRequest);
```

Visual Basic Sample Code

```
If SpectraPro_COM.setMonoSlitWidth(curmono, 1, newWidth) = 0 Then  
    MsgBox ("Error : Failed to Set Slit 1 Width")  
End if
```

Delphi Sample Code

```
if SpectraProServer.setMonoSlitWidth(curmono,1,newWidth) = 0  
  then ShowMessage('Error : Failed to Set Slit 1 Width');
```


MonoSetupScan

Description: Setup the parameters for scanning a monochromator at a fixed rate.

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the connection between the monochromator and NCL (1 or 2).

StartWave: Wavelength to start the scan at

StopWave: Wavelength to end the scan at

ScanRate: Rate in nm per minute to scan at

NumScans: Number of times the scan is to be repeated

SecDelay: Delay in seconds between repeating scans

Result:
0 An invalid request
!0 A valid request

Definitions

Delphi

```
Function MonoSetupScan( MonoNum: Integer;  
                        StartWave: Double;  
                        StopWave: Double;  
                        ScanRate: Double;  
                        NumScans: Integer;  
                        SecDelay: Integer): Integer; safecall;
```

Microsoft IDL

```
[id(0x0000001f)]  
HRESULT MonoSetupScan(  
    [in] long MonoNum,  
    [in] double StartWave,  
    [in] double StopWave,  
    [in] double ScanRate,  
    [in] long NumScans,  
    [in] long SecDelay,  
    [out, retval] long* ValidRequest);
```

Visual Basic Sample Code

```
Dim newStart As Double
Dim newStop As Double
Dim newRate As Double
Dim newNum As Long
Dim newDelay As Long

' set the values
If SpectraPro_COM.MonoSetupScan(curmono, newStart, newStop, newRate, newNum, newDelay) = 0 Then
    MsgBox ("Error : Failed to set Params")
End If
```

Delphi Sample Code

```
var
    newStart,newStop,newRate : double;
    newNum,newDelay : integer;

// set the values
if SpectraProServer.MonoSetupScan(CurMono,newStart,newStop,newRate,newNum,newDelay) = 0
    then showmessage('Error : Failed to set Params');
```

getMonoScanParams

Description: Retrieve the current scan parameters of the selected monochromator

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the connection between the monochromator and NCL (1 or 2).

StartWave: Wavelength to start the scan at

StopWave: Wavelength to end the scan at

ScanRate: Rate in nm per minute to scan at

NumScans: Number of times the scan is to be repeated

SecDelay: Delay in seconds between repeating scans

Result: 0 An invalid request
!0 A valid request

Definitions

Delphi

```
Function getMonoScanParams(MonoNum: Integer;  
                           out StartWave: Double;  
                           out StopWave: Double;  
                           out ScanRate: Double;  
                           out NumScans: Integer;  
                           out SecDelay: Integer): Integer; safecall;
```

Microsoft IDL

```
[id(0x0000001d)]  
HRESULT getMonoScanParams(  
    [in] long MonoNum,  
    [out] double* StartWave,  
    [out] double* StopWave,  
    [out] double* ScanRate,  
    [out] long* NumScans,  
    [out] long* SecDelay,  
    [out, retval] long* ValidRequest);
```

Visual Basic Sample Code

```
Dim newStart As Double
Dim newStop As Double
Dim newRate As Double
Dim newNum As Long
Dim newDelay As Long

If SpectraPro_COM.getMonoScanParams(curmono,newStart,newStop,newRate,newNum,newDelay) <> 0 Then
    ScanStartEdt.Text = newStart
    ScanStopEdt.Text = newStop
    ScanRateEdt.Text = newRate
    ScanNumEdt.Text = newNum
    ScanDelayEdt.Text = newDelay
End if
```

Delphi Sample Code

```
var
newStart,newStop,newRate : double;
newNum,newDelay : integer;

if SpectraProServer.getMonoScanParams(curMono,newStart,newStop,newRate,newNum,newDelay) <> 0
then begin
    ScanStartEdt.Text := format('%.3f',[newStart]);
    ScanStopEdt.Text := format('%.3f',[newStop]);
    ScanRateEdt.Text := format('%.3f',[newRate]);
    ScanNumEdt.Text := IntToStr(newNum);
    ScanDelayEdt.Text := IntToStr(newDelay);
End;
```

MonoScanWavelength

Description: Scan a monochromator between two wavelengths, based on *MonoSetupScan* Values

Parameters

Mono_Num: The number of the monochromator to which commands are addressed.
Note: If you are not using an NCL then use monochromator 1. When using an NCL use the number of the connection between the monochromator and NCL (1 or 2).

Result: 0 An invalid request
!0 A valid request

Definitions

Delphi

Function MonoScanWavelength(MonoNum: Integer): Integer; safecall;

Microsoft IDL

```
[id(0x00000020)]
HRESULT MonoScanWavelength(
    [in] long MonoNum,
    [out, retval] long* ValidResult);
```

Note 2: MonoSetupScan should be called at least once before invoking this function

Visual Basic Sample Code

```
If SpectraPro_COM.MonoScanWavelength(curmono) = 0 Then
    MsgBox ("Error : Failed to Scan")
End If
```

Delphi Sample Code

```
if SpectraProServer.MonoScanWavelength(curMono) = 0
    then showMessage('Error : Failed to Scan');
```

NCL_Present

Description: Determines if an NCL is in the system

Parameters

Result:

0	No NCL present
!0	NCL present

Definitions

Delphi

Function NCL_Present: Integer; safecall;

Microsoft IDL

```
[id(0x00000009)]  
HRESULT NCL_Present([out, retval] long* NCLState);
```

Visual Basic Sample Code

```
If SpectraPro_COM.NCL_Present = 0 Then  
    Labell.Caption = "NCL is not Present"  
End if
```

Delphi Sample Code

```
// display some NCL Info  
if SpectraProServer.NCL_Present = 0  
    then Labell.caption := 'NCL is not Present';
```

Filter_Present

Description: Determines if a filter wheel is defined in the hardware configuration

Parameters

Filter_Num: On systems without an NCL only a value of 1 is valid. On NCL based systems 1 or 2 are valid values, however 2 is reserved for special applications

Result:

0	An invalid request
!0	A valid request

Definitions

Delphi function Filter_Present(Filter_Num: Integer): Integer; safecall;

Microsoft IDL

```
[id(0x0000000b)]
HRESULT Filter_Present(
    [in] long FilterNum,
    [out, retval] long* FilterState);
```

Visual Basic Sample Code

```
If SpectraPro_COM.Filter_Present <> 0 Then
    Labell.caption = "Filter Present"
End If
```

Delphi Sample Code

```
if SpectraProServer.Filter_Present <> 0
then Labell.caption := 'Filter Present';
```

getFilterPosition

Description: Returns the current filter position number in the light path

Parameters

Filter_Num: On systems without an NCL only a value of 1 is valid. On NCL based systems 1 or 2 are valid values, however 2 is reserved for special applications

Result:
0 An invalid request
1- 6 The filter wheel position that is in the light path

Definitions

Delphi function getFilterPosition(Filter_Num: Integer): Integer;
safecall;

:

Microsoft IDL

```
[id(0x00000018)]  
HRESULT getFilterPosition(  
    [in] long FilterNum,  
    [out, retval] long* curFilter);
```

Visual Basic Sample Code

```
If SpectraPro_COM.getFilterPosition(1) = 6 Then  
    Labell.caption = "Filter number 6 is in place"  
End if
```

Delphi Sample Code

```
if SpectraProServer.getFilterPosition(1) = 2  
then Labell.caption := 'filter number 2';
```


setFilterPositon

Description: Sets the filter wheel to a specific position. An invalid position is ignored.

Parameters

Filter_Num: On systems without an NCL only a value of 1 is valid. On NCL based systems 1 or 2 are valid values, however 2 is reserved for special applications

newPosition: Sets the filter wheel to this position

Result: 0 An invalid request
!0 A valid request

Definitions

Delphi function setFilterPosition(Filter_Num: Integer;
NewPosition: Integer
): Integer; safecall;

:

Microsoft IDL

```
[id(0x00000019)]  
HRESULT setFilterPosition(  
    [in] long FilterNum,  
    [in] long newPosition,  
    [out, retval] long* ValidRequest);
```

Visual Basic Sample Code

```
If SpectraPro_COM.setFilterPosition(1, newFilterPos) = 0 Then  
    MsgBox ("Error : Failed to set filter position")  
End If
```

Dephi Sample Code

```
if SpectraProServer.setFilterPosition(1,newFilterPos) = 0  
    then showmessage('Error : Failed to set filter position');
```

FilterHome

Description: Centers the filter on the home (Filter 1) position.

Parameters

Filter_Num: On systems without an NCL only a value of 1 is valid. On NCL based systems 1 or 2 are valid values, however 2 is reserved for special applications

Result: 0 An invalid request
!0 A valid request

Definitions

Delphi function FilterHome(Filter_Num: Integer): Integer; safecall;

Microsoft IDL

```
[id(0x0000001a)]
HRESULT FilterHome(
    [in] long FilterNum,
    [out, retval] long* ValidRequest);
```

Visual Basic Sample Code

```
If SpectraPro_COM.Filter_Present <> 0 Then
    MsgBox ("Homed the Filter wheel")
End If
```

Delphi Sample Code

```
if SpectraProServer.Filter_Present <> 0
    then showmessage('Homed the Filter wheel');
```