

Microbolometer Setup Instructions

1. FPGA cable assignments are located with the cables.
2. The microbolometer is a Photon 320.
3. The VHDL project to use is 320_ethernet
4. The FPGA board to use is the ML-506.
5. The bitstream file to use is already present in the 320_ethernet directory. The VHDL does not need to be recompiled; the FPGA merely needs to be programmed. Programming instructions are included with the project documentation.
6. The Ethernet driver does not work.
7. A Camera Link cable is used to connect the FPGA to the computer framegrabber. A breakout cable is located between the FPGA and the Camera Link cable.
8. The Camera Link output is 324 columns x 256 rows x 14 bits, monochrome.
9. For the Matlab display and acquisition software in the project, use Matlab R2010b SP2, 32 bit version.
10. The camera is connected to the cable labeled "Photon 320 camera". The serial port on the camera cable is NOT to be connected to the serial port on the FPGA board.
11. The configuration dip switches on the FPGA board must be set for the 320. These settings are in the project documentation.
12. The FPGA is connected to the Camera Link cable using the cable labeled "Photon 320 Camera Link"

FPGA Wiring Instructions

1. Connector headers go in the differential column of J4 on the ML-506 board.
2. The header from the Photon 320 breakout cable is plugged in such that the wires are connected to the differential column of J4, with the port marked "22" connected to pin 2 and the other ports connected in increasing order.

Camera Link Breakout Cable

The connector headers on the Camera Link breakout cable are connected to the differential column of J4, as shown in the chart below.

Twisted Pair Colors	Wire 1/Pin 1	Wire 2/Pin 2	Notes
Red, White	White/18	Red/20	Arrow at position 20
Green, Red	Red/22	Green/24	Arrow at position 24
Yellow, Brown	Brown/26	Yellow/28	Arrow at position 28
Pink, Black	Black/30	Pink/32	Arrow at position 32
Brown, White	White/34	Brown/36	Arrow at position 36

To use the chart, chose a wire pair listed in column 1. Plug in the connector such that the arrow is at the position located in column 4. The individual wires should thereby be connected as shown in columns 2 and 3.





3M™ Mini D Ribbon (MDR) Connectors

.050" IDC Wiremount Panelmount Receptacle - Shielded


102 Series



Contact Quantity	Dimensions						
	A	B	C	D	E	F	G
20	.450 [11.43]	1.081 [27.45]	1.31 [33.4]	.65 [16.4]	.92 [23.4]	.71 [18.0]	.94 [23.8]
26	.600 [15.24]	1.231 [31.26]	1.46 [37.2]	.80 [20.2]	1.07 [27.2]	.86 [21.8]	1.09 [27.6]
36	.850 [21.59]	1.481 [37.61]	1.71 [43.5]	1.05 [26.6]	1.32 [33.5]	1.11 [28.1]	1.33 [33.9]
40	.950 [24.13]	1.581 [40.15]	1.81 [46.1]	1.15 [29.1]	1.42 [36.1]	1.21 [30.7]	1.44 [36.5]
50	1.200 [30.48]	1.831 [46.50]	2.06 [52.4]	1.40 [35.5]	1.67 [42.4]	1.46 [37.0]	1.69 [42.8]
68	1.650 [41.91]	2.281 [57.93]	2.51 [63.8]	1.85 [46.9]	2.12 [53.8]	1.91 [48.4]	2.13 [54.2]



Recommended Panel Cut-Out



Note: Panel Thickness .079 [2.00] max.

Inch		[mm]	
Tolerance Unless Noted			
	.0	.00	.000
Inch	±.1	±.01	±.005

() Dimensions for Reference Only

Ordering Information

102XX-0210EC

Contact Quantity
(See Table)

Application tools: 10951 Assembly Press and 10952 Cable Cutter.

TS-0425-12
Sheet 2 of 2



Interconnect Solutions
<http://www.3M.com/interconnects/>

3M is a trademark of 3M Company.
For technical, sales or ordering information call
800-225-5373

Miscellaneous, Scratch, and Superseded Documentation

Microblaze Setup Instructions

P1
14/Nov 2013

- 1, Cable setups are on sheets with the cable.
- 2, The camera is a Photon 320
- 3, The VHDL project to use is 320_ethernet
- 4, The FPGA board to use is the ML-506.
- 5, The bitstream file to use is already present in the 320_ethernet directory. Programming instructions are with the project documentation
- 6, The ethernet driver does not work.
- 7, A ~~ca~~ Camera Link cable is used to connect the FPGA to the computer framegrabber.
- 8, the output is 324 columns x 256 rows x 14 bits
- 9, Use Matlab[®] R2010b SP2 32 bit version
- 10, Connect to the camera the cable labeled "Photon 320 camera"
- 11, ~~ca~~ ^{ca} connect the FPGA to the ~~ca~~ Camera Link cable using the ~~cable~~ cable labeled "Photon 320 Camera Link"

Photon 320 CameraLink Cable Pinout

pl
12 Nov 2013

~~Break~~
Connector headers go in the differential column of J4 on the ML506 board

Twisted Pair Colors	Wire 1/Pin 1	Wire 2/Pin 2	Notes
Red, White	White/18	Red/20	Arrow at pos 20
Green, Red	Red/22	Green/24	Arrow at pos 24
Yellow, Brown	Brown/26	Yellow/28	Arrow at pos 28
Pink, Black	Black/30	Pink/32	Arrow at pos 32
Brown, White	White/34	Brown/36	Arrow at pos 36

To use, choose a wire pair listed in column 1. Plug in the connector such that the arrow is at the position located in column 4. The individual wires should then be connected as shown in columns 2 and 3.

Addendum Photon 320 Breakout Cable

Plug in the ~~break~~ Photon 320 breakout cable such that the wires are connected to the differential column of J4, with the part marked "21" connected to pin 2 and the other pins in increasing order.

Connection	Wire #	Header Group	Header Pin	FPGA Pin (ML402)	FPGA Pin (ML506)
Synchronization -	12	J4 on ML506	22	Y18	T34
Synchronization +	2		24	AA18	U33
CLK -	11		26	W19	U31
CLK +	7		28	Y19	U32
Data -	13		30	Y21	V33
Data +	8		32	Y20	V32

Header Pin	Header Pin	Net Name
HDR2-22-SM11-N	22	rx Data 0V
HDR2-24-SM11-P	24	rx Data 0P
HDR2-26-SM11-N	26	rx Data 0N
HDR2-28-SM11-P	28	rx Data 0P
HDR2-30-DIFF-3V	30	rx Data 0P
HDR2-32-DIFF-3P	32	rx Data 0P

Pin 6 analog vid ~~rx~~ -
 Pin 1 analog vid ~~rx~~ +

Gnd 3, 4, 5,
 Power 1.5, 10

TX 9 - BDB 9 Female 2
 RX 14 - DB 9 3

~~DTE TX~~
 2p DTE TX Pin 3
 RX Pin 2

P32	HDR2-12-DIFF-0P	12	+data
N32	HDR2-10-DIFF-0-N	10	IBdat
L34	HDR2-4-SM11-2-P	4	ib
K34	HDR2-2-SM11-4	2	ips
K33	HDR2-8-SM11-P	8	
K32	HDR2-6-SM11-7-P	6	ib

Header Pin	Net Name
2	ips
4	ib
6	ib
8	i
10	Idat
12	IBdat


FPGA is DCE wired

Docs say null modem cable needed for computer serial link, but comps are generally DTE, not DCE
 Schematic says DTE, and has pin 2 hooked up to RX, and pin labeled as in.

FIFO = <x> means x is written to the FIFO


data
MOSFET
Suppressed

WREN is a Mealy output. Counter is incremented in a Mealy way



clk-	5	R32	ID-LEAF-SAMPLE II	Reader Pin	Header Name
clk+	18	R32	INTERNAL-SAMPLE II	18	ADDR1-18-DIFF-2M
X0-	2	R34	INTERNAL-SAMPLE II	20	ADDR1-20-DIFF-1P
X0+	15	U32	INTERNAL-SAMPLE II	22	ADDR1-22-SM-L-N
X1-	3	U31	INTERNAL-SAMPLE II	24	ADDR1-24-SM-L-P
X1+	16	U32	INTERNAL-SAMPLE II	26	ADDR1-26-SM-L-N
X2-	4	U32, U33	INTERNAL-SAMPLE II	28	ADDR1-28-SM-L-P
X2+	17	U32, U32	INTERNAL-SAMPLE II	30	
X3-	6	U34	INTERNAL-SAMPLE II	32	
X3+	19	U34	INTERNAL-SAMPLE II	34	
			INTERNAL-SAMPLE II	36	

Change output outputs on state transition from C to A



Entity
 cam-clk in
 cam-data in
 start in
 resetn in
 data-valid buffer
 frame valid buffer
 write-out
 write-en out
 dirty-out
 finished out
 reset-out out
 deser-data buffer state register (IS DOWN TO 0)