LEO 906
Pumping System
Model SQ 095

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BEDIENUNGSHANDBUCH

INSTRUCTION MANUAL

VARIAN

vacuum technologies

Dear Customer,

Thank you for purchasing a VARIAN vacuum product. At VARIAN Vacuum Technologies we make every effort to ensure that you will be satisfied with the product and/or service you have purchased.

As part of our Continuous Improvement effort, we ask that you report to us any problem you may have had with the purchase or operation of our product. On the back side you find a Corrective Action Request form that you may fill out in the first part and return to us.

This form is intended to supplement normal lines of communications and to resolve problems that existing systems are not addressing in an adequate or timely manner.

Upon receipt of your Corrective Action Request we will determine the Root Cause of the problem and take the necessary actions to eliminate it. You will be contacted by one of our employees who will review the problem with you and update you, with the second part of the same form, on our actions.

Your business is very important to us. Please, take the time and let us know how we can improve.

Sincerely,

Sergio PIRAS

*Vice President and General Manager
VARIAN Vacuum Technologies*

Note: Fax or mail the Customer Request for Action (see backside page) to VARIAN Vacuum Technologies (Torino) - Quality Assurance or to your nearest VARIAN representative for onward transmission to the same address.

CUSTOMER REQUEST FOR CORRECTIVE / PREVENTIVE / IMPROVEMENT ACTION

TO : VARIAN VACUUM TECHNOLOGIES TORINO - QUALITY ASSURANCE

FAX N° : XXXX - 011 - 9979350

ADDRESS: VARIAN S.p.A. - Via F.lli Varian, 54 - 10040 Leinì (Torino) - Italy

E-MAIL : marco.marzio@varianinc.com

NAME _____	COMPANY _____	FUNCTION _____
<p>ADDRESS : _____</p> <p>TEL. N° : _____ FAX N° : _____</p> <p>E-MAIL : _____</p>		
<p>PROBLEM / SUGGESTION :</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>		
<p>REFERENCE INFORMATION (model n°, serial n°, ordering information, time to failure after installation, etc.) :</p> <p>_____</p> <p>_____</p> <p>_____</p> <p style="text-align: right;">DATE _____</p>		

<p>CORRECTIVE ACTION PLAN / ACTUATION (by VARIAN VTT)</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>LOG N° _____</p>
--	---------------------

XXXX = Code for dialing Italy from your country (es. 01139 from USA; 00139 from Japan, etc.)

ISTRUZIONI PER L'USO.....	1
GEBRAUCHSANLEITUNG	5
INSTRUCTIONS FOR USE.....	9

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
THIS IS A COMPUTER GENERATED DRAWING

CHANGED FROM		CLASS	
S	F	DIMENSION	H
F	F	MACHINING	
Q	Q	WELD CONSTRUCTION	
C46-94		TV 250 BASE PCB MAINS WIR DIAG	
NUMBER	POSITION	DRAWN	CHECKED
C.M.	C.M.	C.M.	C.M.
DATE	DATE	SCALE	SCALE
06/94	06/94	A3	A
REVISION	REVISION	DRAWING No.	REVISION
A	A	03.660941	A

NOTES

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GENERAL INFORMATION

This equipment is destined for use by professionals. The user should read this instruction manual and any other additional information supplied by Varian before operating the equipment. Varian will not be held responsible for any events occurring due to non-compliance, even partial, with these instructions, improper use by untrained persons, non-authorized interference with the equipment or any action contrary to that provided for by specific national standards. The following paragraphs contain all the information necessary to guarantee the safety of the operator when using the equipment. Detailed information about the installed equipment are available into the supplied relevant technical manuals.

This manual uses the following standard protocol:

! WARNING!

The warning messages are for attracting the attention of the operator to a particular procedure or practice which, if not followed correctly, could lead to serious injury.

☞ CAUTION

The caution messages are displayed before procedures which, if not followed, could cause damage to the equipment.

NOTE



The notes contain important information taken from the text.

DESCRIPTION

The LEO 906 is a pumping group and is composed of the following units:

- an assembly that comprises the Turbo-V 300 SF pump, two vibration dampers and the vacuum chamber with a Cold Cathode Gauge installed;
- an accessories set.

The following figures show the pump and vacuum chamber assembly view, and an overall view with all the accessories mounted.

The following table, with reference to the previous figures, lists the LEO 906 components (N.I. = not indicated).

MOD.	POS.	ITEM
969-9068	1	Turbo-V300SF pump
03.662053	2	Vacuum chamber
03.662200	3	Vibration dampers
03.662315	4	Plate
969-9833	5	Vent valve
27.179700-01	6	KF16 24V Valve
27.179706-01	7	KF25 24V Valve
27.179710-01	8	KF40 24V Valve
27.179705-01	9	KF25 24V Valve
27.229975-01	10	Cold Cathode Gauge
27.229976-01	11	Thermovac Gauge
03.662218	12	KF25/KF16 collector
03.662220	13	KF40/KF25 collector
AKF08-40.25SS	14	KF40/25 adapter
AKF08-25.16SS	15	KF25/16 adapter
AKF14-25.10SS	16	KF25 LG 100 extension
ZTF-40.035SS	17	Metal flexible KF40 LG 350
AKF10-25.050SS	18	Metal flexible KF25 LG 50
AKF10-25.050SS	19	Metal flexible KF25 LG 1000
AKF10-16.050SS	20	Metal flexible KF16 LG 500
AKF10-25.025SS	21	Metal flexible KF25 LG 250
AKF01-32.40AL	22	Clamp KF32/40
AKF01-20.25AL	23	Clamp KF20/25
AKF01-10.16AL	24	Clamp KF10/16
969-9425	N.I.	Controller
03.662054	N.I.	Cryopanel

Remote start/stop, pump status signals, forepump start/stop, interlock control (for pressure switch, water flow switch, etc.) capability, are provided via the controller auxiliary connectors.


STORAGE

When transporting and storing the LEO 906, the following environmental requirements should be satisfied:

- temperature: from -20 °C to +70 °C
- relative humidity: 0 - 95% (without condensation)

PREPARATION FOR INSTALLATION

The LEO 906 is supplied in two special protective packings: the first one contains the assembly pump-column, the second one contains all the other accessories. If these show signs of damage which may have occurred during transport, contact your local sales office. When unpacking the LEO 906, ensure that the system is not dropped or subjected to any form of impact.

Pump-column assembly packaging

Do not dispose of the packing materials in an unauthorised manner. The material is 100% recyclable and complies with EEC Directive 85/399.

INSTALLATION

WARNING!

Cause its weight, the LEO 906 must be handled by means of suitable moving and handling tools.

The LEO 906 installation must be executed according to what stated in the previous figures.

The pump-column assembly is put into a shipping structure that must be removed by the customer before the installation.

To remove the structure four M& screws must be unscrewed: two screws fix the assembly base plate to two structure brackets, two screws fix the structure to two brackets welded to the column.

Four M6 screws must be used in the final installation to fix the assembly to suitable brackets.

WARNING!

Turbo Stack must be fixed to the standing plane by means of suitable fixing device.

Detailed information about the other units installed into the LEO 906 are contained in the relevant Instructions Manuals.

WARNING!

The Turbo pump of the LEO 906 can reach high temperature that can cause severe injury. Pay attention during the unit handling.

When the Turbo pump is water cooled, be sure that the water tubes are not in contact with live elements.

NOTE

The LEO 906 installed in its final position must be positioned so that cold air can flow through the units. Do not install or use the LEO 906 in an environment exposed to atmospheric agents (rain, snow, ice), dust, aggressive gases, or in explosive environments or those with a high fire risk.

During operation, the following environmental conditions must be respected:

- temperature: from 0 °C to +40 °C
- relative humidity: 0 - 95% (without condensation)

USE

WARNING!

When employing the LEO 906 for pumping toxic, flammable, or radioactive gases, please follow the required procedures for each gas disposal.

WARNING!

Never operate the LEO 906 if the inlet flanges are not connected to the system or blanked off.

WARNING!

When the pump must be reinstalled after a removal, to avoid dangers in the case of a sudden stop of the Turbo pump, it must be fixed to the LEO 906 by means of the suitable fixing device using screws with a minimum yield point of 500 N/mm² and fixing them with a torque of 3.1 Nm.

Detailed information and operating procedures that involve optional connections or options are supplied in the Instruction Manual of each unit.

MAINTENANCE

The LEO 906 does not require any maintenance. Any work performed on the LEO 906 must be carried out by authorised personnel.

When an unit mounted on the LEO 906 is faulty, refer to the relevant Instruction Manual.

WARNING!


Before carrying out any work on the LEO 906, disconnect it from the supply.

WARNING!

Should the LEO 906 be used for pumping toxic, flammable, or radioactive gases and it is to be scrapped, please follow the required procedures for each gas disposal.

The following figure shows the LEO 906 block diagram.

THIS IS A COMPUTER GENERATED DRAWING

CHANGED FROM		CLASS	
S.F.	F.O.	H	
C46-94	NUMBER	TOLERANCES (Unless otherwise specified)	
C.M.	POSITION	DIMENSION	±
06/94	DATE	X0 X30 X720 X315 X1000X2000	ANGLE ±
A	REVISION	X6 X20 X75 X100 X200	PERPENDICULARITY
		MACHINING ±0.2 ±0.3 ±0.5 ±0.8 ±1.2 ±2	PARALLELISM
		WELD CONSTRUCTION ±0.5 ±0.8 ±1 ±2 ±3 ±4	CHAMFER
		OTHER	
TV 250 BASE PCB MAINS WIR DIAG			
JOB ORDER		DRAWN	CHECKED
C.M.		C.M.	C.M.
DATE		DATE	SCALE
06/94		07/06/94	
REVISED		SIZE	REVISION
A		A3	A
		DRAWING No.	
		03.660941	

NOTES

Block diagram

Request for Return

1. A Return Authorization Number (RA#) **WILL NOT** be issued until this Request for Return is completely filled out, signed and returned to Varian Customer Service.
2. Return shipments shall be made in compliance with local and international **Shipping Regulations** (IATA, DOT, UN).
3. The customer is expected to take the following actions to ensure the **Safety** of workers at Varian: (a) Drain any oils or other liquids, (b) Purge or flush all gasses, (c) Wipe off any excess residues in or on the equipment, (d) Package the equipment to prevent shipping damage, (for Advance Exchanges please use packing material from replacement unit).
4. Make sure the shipping documents clearly show the RA# and then return the package to the Varian location nearest you.

North and South America

Varian Vacuum Technologies
 121 Hartwell Ave
 Lexington, MA 02421
 Phone : +1 781 8617200
 Fax: +1 781 8609252

Europe and Middle East

Varian SpA
 Via Flli Varian 54
 10040 Leini (TO) – ITALY
 Phone: +39 011 9979111
 Fax: +39 011 9979330

Asia and ROW

Varian Vacuum Technologies
 Local Office

CUSTOMER INFORMATION

Company name:	
Contact person: Name:	Tel:
Fax:	E-Mail:
Ship Method:	Shipping Collect #: P.O.#:
<i>Europe only:</i> VAT reg. Number:	<i>USA only:</i> <input type="checkbox"/> Taxable <input type="checkbox"/> Non-taxable
Customer Ship To:	Customer Bill To:
.....
.....

PRODUCT IDENTIFICATION

Product Description	Varian P/N	Varian S/N	Purchase Reference

TYPE OF RETURN (check appropriate box)

<input type="checkbox"/> Paid Exchange	<input type="checkbox"/> Paid Repair	<input type="checkbox"/> Warranty Exchange	<input type="checkbox"/> Warranty Repair	<input type="checkbox"/> Loaner Return
<input type="checkbox"/> Credit	<input type="checkbox"/> Shipping Error	<input type="checkbox"/> Evaluation Return	<input type="checkbox"/> Calibration	<input type="checkbox"/> Other

HEALTH and SAFETY CERTIFICATION

Varian Vacuum Technologies **CAN NOT ACCEPT** any equipment which contains **BIOLOGICAL HAZARDS** or **RADIOACTIVITY**. Call Varian Customer Service to discuss alternatives if this requirement presents a problem.

The equipment listed above (check one):

HAS NOT been exposed to any toxic or hazardous materials

OR

HAS been exposed to any toxic or hazardous materials. In case of this selection, check boxes for any materials that equipment was exposed to, check all categories that apply:

Toxic Corrosive Reactive Flammable Explosive Biological Radioactive

List all toxic or hazardous materials. Include product name, chemical name and chemical symbol or formula.

.....

Print Name: Customer Authorized Signature:

Print Title: Date:/...../.....

NOTE: If a product is received at Varian which is contaminated with a toxic or hazardous material that was not disclosed, **the customer will be held responsible** for all costs incurred to ensure the safe handling of the product, and **is liable** for any harm or injury to Varian employees as well as to any third party occurring as a result of exposure to toxic or hazardous materials present in the product.

Do not write below this line

Notification (RA)#: Customer ID#: Equipment #:

FAILURE REPORT

TURBO PUMPS and TURBOCONTROLLERS

<input type="checkbox"/> Does not start <input type="checkbox"/> Does not spin freely <input type="checkbox"/> Does not reach full speed <input type="checkbox"/> Mechanical Contact <input type="checkbox"/> Cooling defective	<input type="checkbox"/> Noise <input type="checkbox"/> Vibrations <input type="checkbox"/> Leak <input type="checkbox"/> Overtemperature	POSITION <input type="checkbox"/> Vertical <input type="checkbox"/> Horizontal <input type="checkbox"/> Upside-down <input type="checkbox"/> Other:	PARAMETERS Power: Rotational Speed: Current: Inlet Pressure: Temp 1: Foreline Pressure: Temp 2: Purge flow: <hr/> OPERATION TIME:
TURBOCONTROLLER ERROR MESSAGE:			

ION PUMPS/CONTROLLERS

<input type="checkbox"/> Bad feedthrough <input type="checkbox"/> Vacuum leak <input type="checkbox"/> Error code on display	<input type="checkbox"/> Poor vacuum <input type="checkbox"/> High voltage problem <input type="checkbox"/> Other
Customer application:	

VALVES/COMPONENTS

<input type="checkbox"/> Main seal leak <input type="checkbox"/> Solenoid failure <input type="checkbox"/> Damaged sealing area	<input type="checkbox"/> Bellows leak <input type="checkbox"/> Damaged flange <input type="checkbox"/> Other
Customer application:	

LEAK DETECTORS

<input type="checkbox"/> Cannot calibrate <input type="checkbox"/> Vacuum system unstable <input type="checkbox"/> Failed to start	<input type="checkbox"/> No zero/high background <input type="checkbox"/> Cannot reach test mode <input type="checkbox"/> Other
Customer application:	

INSTRUMENTS

<input type="checkbox"/> Gauge tube not working <input type="checkbox"/> Communication failure <input type="checkbox"/> Error code on display	<input type="checkbox"/> Display problem <input type="checkbox"/> Degas not working <input type="checkbox"/> Other
Customer application:	

PRIMARY PUMPS

<input type="checkbox"/> Pump doesn't start <input type="checkbox"/> Doesn't reach vacuum <input type="checkbox"/> Pump seized	<input type="checkbox"/> Noisy pump (describe) <input type="checkbox"/> Over temperature <input type="checkbox"/> Other
Customer application:	

DIFFUSION PUMPS

<input type="checkbox"/> Heater failure <input type="checkbox"/> Doesn't reach vacuum <input type="checkbox"/> Vacuum leak	<input type="checkbox"/> Electrical problem <input type="checkbox"/> Cooling coil damage <input type="checkbox"/> Other
Customer application:	

FAILURE DESCRIPTION

(Please describe in detail the nature of the malfunction to assist us in performing failure analysis):

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India

Varian India PVT LTD

101-108, 1st Floor
1010 Competent House
7, Nangal Raya Business Centre
New Delhi 110 046
India
Tel: (91) 11 5548444
Fax: (91) 11 5548445

Italy

Varian Vacuum Technologies

Via F.lli Varian, 54
10040 Leini, (Torino)
Italy
Tel: (39) 011 997 9111
Fax: (39) 011 997 9350

Japan

Varian Vacuum Technologies

Sumitomo Shibaura Building, 8th Floor
4-16-36 Shibaura
Minato-ku, Tokyo 108
Japan
Tel: (81) 3 5232 1253
Fax: (81) 3 5232 1263

Korea

Varian Technologies Korea, Ltd.

Shinsa 2nd Bldg. 2F
966-5 Daechi-dong
Kangnam-gu, Seoul
Korea 135-280
Tel: (82) 2 3452 2452
Fax: (82) 2 3452 2451

Mexico

Varian S.A.

Concepcion Beistegui No 109
Col Del Valle
C.P. 03100
Mexico, D.F.
Tel: (52) 5 523 9465
Fax: (52) 5 523 9472

Taiwan

Varian Technologies Asia Ltd.

18F-13 No.79, Hsin Tai Wu Road
Sec. 1, Hsi Chih
Taipei Hsien
Taiwan, R.O.C.
Tel: (886) 2 2698 9555
Fax: (886) 2 2698 9678

UK and Ireland

Varian Ltd.

28 Manor Road
Walton-On-Thames
Surrey KT 12 2QF
England
Tel: (44) 1932 89 8000
Fax: (44) 1932 22 8769

United States

Varian Vacuum Technologies

121 Hartwell Avenue
Lexington, MA 02421
USA
Tel: (781) 861 7200
Fax: (781) 860 5437
Toll Free: (800) 882 7426

Other Countries

Varian Vacuum Technologies

Via F.lli Varian, 54
10040 Leini, (Torino)
Italy
Tel: (39) 011 997 9111
Fax: (39) 011 997 9350

Internet Users:

Customer Service & Technical Support:

vtt.customer.service@varianinc.com

Worldwide Web Site:

www.varianinc.com/vacuum

Order On-line:

www.evarian.com

Representatives in most countries

VARIAN